

StimConnect[®]

Disposable connection cables for PACE[®] Model 203 (P/N 2002)

StimConnect[®] connection cables are available in color and length as indicated.



PACE[®] Model 203 with collets (clamping terminals) accepting pins of 0.9 to 2.0 mm

(order number: P/N 2002)







Connection pacemaker	Disposable connection cable StimConnect [®]		Connection pacing electrode/ lead
	Connection to pacemaker	Connection to pacing electrode	
	Osypka plug with 2 mm pins, touch protected Length: 2.0 m SC.B/TSP/OSY (blue)	Terminal for shrouded 2 mm pins Order number: SC.W/TSP/OSY (white)	Example: Plugs with unipolar shrouded pin
	Osypka plug with 2 mm pins, touch protected Length: 2.0 m SC.B/TBOX/OSY (blue)	Terminal with sided clamping screws for exposed 0.1 to 2.0 mm pins Order number: SC.W/TBOX/OSY (white)	Example: Plugs with unipolar exposed pin
	Osypka plug with 2 mm pins, touch protected Length: 2.5 m SC.B/TCO/OSY (blue)	Terminal for exposed 0.9 to 2.0 mm pins Order number: SC.W/TCO/OSY (white)	Example: Plugs with unipolar exposed pin
	Plug with unipolar exposed 2 mm pins, with protection cap Length: 2.0 m SC.B/TSP/OP (blue)	Terminal for shrouded 2 mm pins Order number: SC.W/TSP/OP (white)	Example: Plugs with unipolar shrouded pin
	Plug with unipolar exposed 2 mm pins, with protection cap Length: 2.0 m SC.B/TBOX/OP (blue)	Terminal with sided clamping screws for exposed 0.1 to 2.0 mm pins Order number: SC.W/TBOX/OP (white)	Example: Plugs with unipolar exposed pin
	Plug with unipolar exposed 2 mm pins, with protection cap Length: 2.5 m SC.W/TCO/OP (white)	Terminal for exposed 0.9 to 2.0 mm pins Order number: SC.W/TCO/OP (white)	Example: Plugs with unipolar exposed pin

StimConnect®

Reusable connection cables for PACE® Model 203 (P/N 2002)

StimConnect® connection cables are available in color and length as indicated

Connection pacemaker	Reusable connection cable StimConnect®		Connection pacing electrode/ lead
	Connection to pacemaker	Connection to pacing electrode	
 <p>PACE® Model 203 with collets (clamping terminals) accepting pins of 0.9 to 2.0 mm</p> <p>(order number: P/N 2002)</p>	<p>Plug with unipolar exposed 2 mm pins, with protection cap</p>  <p>Length: 2.0 m Order number: SC.RB/TBOX/SP (blue)</p>	<p>Terminal with sided clamping screws for exposed 0.1 to 2.0 mm pins</p>  <p>Order number: SC.RW/TBOX/SP (white)</p>	<p>Example: Plugs with unipolar exposed pin</p> 

SCPAC203EN220412C