


EndoStim™

TEMPORARY PACING LEADS



Temporary Pacing Leads

EndoStim™

Pacing leads for transvenous applications to temporarily stimulate the ventricle or the atrium. Based on the application of external pacemakers suitable for emergency, pre- and postoperative stimulation as well as for rapid-pacing during TAVI procedures.

Characteristics


The EndoStim™ portfolio includes bipolar and quadripolar pacing leads. A floating balloon catheter completes the product range.

The following French sizes, shapes, lengths, inter-electrode distances and connectors are available:

	Bipolar EndoStim™ Leads	Quadripolar EndoStim™ Leads	Pacing Catheter with Balloon
French Size	4F, 5F, 6F	4F, 5F, 6F	5F
Length	90 cm, 125 cm	90 cm, 125 cm	110 cm
Shape	Straight, Cournand, Josephson	Straight, Josephson	Cournand
Inter-Electrode Distance	2 mm, 5 mm, 10 mm	5 mm, 10 mm, 2-5-2 mm	10 mm
Connector	Exposed 2 mm pin with protective cap	Shrouded 2 mm pin	Exposed 2 mm pin with protective cap

Shape

EndoStim™ pacing leads are available in three shapes: Straight, Cournand and Josephson.



Connector

EndoStim™ pacing leads are either available with shrouded pins or with exposed pins with protective cap.



shrouded pins
(Bipolar pacing leads)



exposed pins with protective cap
(Quadripolar pacing leads, pacing lead with balloon)

Order Information Temporary Pacing Leads

Bipolar EndoStim™ Temporary Pacing Leads				(Packaging Unit: 5 pieces)		
Article No.	Number of Poles	French Size	Length	Shape	Inter-Electrode Distance	Connector
27039-OMG	2	4F	90 cm	Straight	2 mm	Exposed Pins
27003-F	2	4F	90 cm	Straight	5 mm	Exposed Pins
27010-OMG	2	4F	90 cm	Straight	10 mm	Exposed Pins
27040-OMG	2	4F	90 cm	Cournand	2 mm	Exposed Pins
27018-OMG	2	4F	90 cm	Cournand	5 mm	Exposed Pins
27015-OMG	2	4F	90 cm	Cournand	10 mm	Exposed Pins
27038-OMG	2	4F	90 cm	Josephson	2 mm	Exposed Pins
27006-F	2	4F	90 cm	Josephson	5 mm	Exposed Pins
27009-OMG	2	4F	90 cm	Josephson	10 mm	Exposed Pins
27051-OMG	2	5F	90 cm	Straight	2 mm	Exposed Pins
27045-OMG	2	5F	90 cm	Straight	5 mm	Exposed Pins
27048-OMG	2	5F	90 cm	Straight	10 mm	Exposed Pins
27052-OMG	2	5F	90 cm	Cournand	2 mm	Exposed Pins
27046-OMG	2	5F	90 cm	Cournand	5 mm	Exposed Pins
27049-OMG	2	5F	90 cm	Cournand	10 mm	Exposed Pins
27050-OMG	2	5F	90 cm	Josephson	2 mm	Exposed Pins
27044-OMG	2	5F	90 cm	Josephson	5 mm	Exposed Pins
27047-OMG	2	5F	90 cm	Josephson	10 mm	Exposed Pins
27042-OMG	2	5F	125 cm	Straight	2 mm	Exposed Pins
27004-F	2	5F	125 cm	Straight	5 mm	Exposed Pins
27012-OMG	2	5F	125 cm	Straight	10 mm	Exposed Pins
27043-OMG	2	5F	125 cm	Cournand	2 mm	Exposed Pins
27019-OMG	2	5F	125 cm	Cournand	5 mm	Exposed Pins
27016-OMG	2	5F	125 cm	Cournand	10 mm	Exposed Pins
27041-OMG	2	5F	125 cm	Josephson	2 mm	Exposed Pins
27007-F	2	5F	125 cm	Josephson	5 mm	Exposed Pins
27011-OMG	2	5F	125 cm	Josephson	10 mm	Exposed Pins
27063-OMG	2	6F	90 cm	Straight	2 mm	Exposed Pins
27057-OMG	2	6F	90 cm	Straight	5 mm	Exposed Pins
27060-OMG	2	6F	90 cm	Straight	10 mm	Exposed Pins
27064-OMG	2	6F	90 cm	Cournand	2 mm	Exposed Pins
27058-OMG	2	6F	90 cm	Cournand	5 mm	Exposed Pins
27061-OMG	2	6F	90 cm	Cournand	10 mm	Exposed Pins
27062-OMG	2	6F	90 cm	Josephson	2 mm	Exposed Pins
27056-OMG	2	6F	90 cm	Josephson	5 mm	Exposed Pins
27059-OMG	2	6F	90 cm	Josephson	10 mm	Exposed Pins
27054-OMG	2	6F	125 cm	Straight	2 mm	Exposed Pins
27005-F	2	6F	125 cm	Straight	5 mm	Exposed Pins
27014-OMG	2	6F	125 cm	Straight	10 mm	Exposed Pins
27055-OMG	2	6F	125 cm	Cournand	2 mm	Exposed Pins
27020-OMG	2	6F	125 cm	Cournand	5 mm	Exposed Pins
27017-OMG	2	6F	125 cm	Cournand	10 mm	Exposed Pins
27053-OMG	2	6F	125 cm	Josephson	2 mm	Exposed Pins
27008-F	2	6F	125 cm	Josephson	5 mm	Exposed Pins
27013-OMG	2	6F	125 cm	Josephson	10 mm	Exposed Pins

Temporary Pacing Leads

Order Information Temporary Pacing Leads

Quadripolar EndoStim™ Temporary Pacing Leads				(Packaging Unit: 5 pieces)		
Article No.	Number of Poles	French Size	Length	Shape	Inter-Electrode Distance	Connector
27022-OMG	4	4F	90 cm	Straight	5 mm	Shrouded Pins
27028-OMG	4	4F	90 cm	Straight	10 mm	Shrouded Pins
27033-OMG	4	4F	90 cm	Straight	2-5-2 mm	Shrouded Pins
27021-OMG	4	4F	90 cm	Josephson	5 mm	Shrouded Pins
27027-OMG	4	4F	90 cm	Josephson	10 mm	Shrouded Pins
27032-OMG	4	4F	90 cm	Josephson	2-5-2 mm	Shrouded Pins
27024-OMG	4	5F	125 cm	Straight	5 mm	Shrouded Pins
27030-OMG	4	5F	125 cm	Straight	10 mm	Shrouded Pins
27035-OMG	4	5F	125 cm	Straight	2-5-2 mm	Shrouded Pins
27023-OMG	4	5F	125 cm	Josephson	5 mm	Shrouded Pins
27029-OMG	4	5F	125 cm	Josephson	10 mm	Shrouded Pins
27034-OMG	4	5F	125 cm	Josephson	2-5-2 mm	Shrouded Pins
27026-OMG	4	6F	125 cm	Straight	5 mm	Shrouded Pins
27065-OMG	4	6F	125 cm	Straight	10 mm	Shrouded Pins
27037-OMG	4	6F	125 cm	Straight	2-5-2 mm	Shrouded Pins
27025-OMG	4	6F	125 cm	Josephson	5 mm	Shrouded Pins
27031-OMG	4	6F	125 cm	Josephson	10 mm	Shrouded Pins
27036-OMG	4	6F	125 cm	Josephson	2-5-2 mm	Shrouded Pins

Bipolar EndoStim™ Temporary Pacing Lead with Balloon								(Packaging Unit: 1 piece)
Article No.	Number of Poles	French Size	Length	Shape	Inter-Electrode Distance	Ballon Diameter	Ballon Volume	Connector
29207-BS	2	5F	110 cm	Cournand	10 mm	9 mm	1.0 cc	Exposed Pins



EasyStim™ Adapter – Accessory for Temporary Pacing Leads

The EasyStim™ adapter allows to connect shrouded 2mm pins to collets for exposed (naked) 2mm pins.

EasyStim™ Adapter		(Packaging Unit: 20 double packs)
Article No.	Product Name	Sterility
94005	EasyStim	no-sterile
94006	EasyStim.S	sterile



Osypka Medical GmbH

Albert-Einstein-Strasse 3
 12489 Berlin, Germany
 Phone: +49 (30) 6392 8300
 Fax: +49 (30) 6392 8301
 E-Mail: mail@osypkamed.com
 www.osypkamed.com



EasyStim and EndoStim are trademarks of Osypka Medical,
 SPIKE FLOW is a trademark of FIAB
 © Osypka Medical 2022



0344