

# StimConnect™

## for PACE® Model 101 with Unipolar Terminals

### Extension Cables for Transvenous Pacing Applications

PACE® Model 101 is primarily used for temporary transvenous cardiac pacing applications. Transvenous leads for temporary pacing are either equipped with 2 mm shrouded plug pins or exposed (naked) plug pins. Following StimConnect® extension cables are recommended for PACE® Model 101 with unipolar terminals.

Pacemaker	Disposable Extension Cable StimConnect™		Transvenous Pacing Lead
	Connection Pacemaker- Side	Connection Patient-Side	
 <p>PACE® Model 101 with unipolar terminals accepting pins of 0.9 to 2.0 mm (order number: P/N 2001).</p> 	 <p><b>Article No: 82207BL</b></p> <p>Touch protected plug fitting unipolar terminals of PACE® Model 101</p> <p>Connects shrouded pins of transvenous leads</p> <p>Packaging Unit: 10 Pieces</p>	 <p>Shrouded pins (ø 2 mm)</p>	
	 <p><b>Article No: 81987BL</b></p> <p>Touch protected plug fitting unipolar terminals of PACE® Model 101</p> <p>Connects exposed (naked) pins of transvenous leads</p> <p>Packaging Unit: 5 Pieces</p>	 <p>Exposed (naked) pins (ø 2 mm)</p>	

Please contact [mail@osypkamed.com](mailto:mail@osypkamed.com) for other StimConnect® extension cable variants.

# StimConnect™

## for PACE® Model 101 with Bipolar Terminal

### Extension Cables for Transvenous Pacing Applications

PACE® Model 101 is primarily used for temporary transvenous cardiac pacing applications. Transvenous leads for temporary pacing are either equipped with 2 mm shrouded plug pins or exposed (naked) plug pins. Following StimConnect® extension cables are recommended for PACE® Model 101 with bipolar terminals.

Pacemaker	Disposable Extension Cable StimConnect™		Transvenous Pacing Lead
	Connection Pacemaker- Side	Connection Patient-Side	
 <p>PACE® Model 101 with bipolar terminal (order number: P/N 2004)</p>	 <p><b>Article No: 82208BL</b></p> <p>Bipolar plug fitting bipolar terminal of PACE® Model 101</p> <p>Connects shrouded pins of transvenous leads</p> <p>Packaging Unit: 10 Pieces</p>	 <p>Shrouded pins (ø 2 mm)</p>	
	 <p><b>Article No: 82203BL</b></p> <p>Bipolar plug fitting bipolar terminal of PACE® Model 101</p> <p>Connects exposed (naked) pins of transvenous leads</p> <p>Packaging Unit: 10 Pieces</p>	 <p>Exposed (naked) pins (ø 2 mm)</p>	

**OSYPKA**  
MEDICAL







Please contact [mail@osypkamed.com](mailto:mail@osypkamed.com) for other StimConnect® extension cable variants.

# StimConnect™ for PACE® Model 203 with Unipolar Terminals

## Extension Cables for Surgical Pacing Applications

PACE® Model 203 is primarily used for temporary surgical cardiac pacing applications. Surgical wires (temporary myocardial electrodes or heart wires) for temporary pacing are either equipped with 2 mm exposed (naked) plug pins or 0.9 mm pins/ needles. Following StimConnect® extension cables are recommended for PACE® Model 203 with unipolar terminals.

Disposable cables are for single use, reusable extension cables can be sterilized 25 times by ethylene oxide or steam. A punch card allows the identification of the sterilization cycles.

Pacemaker	Disposable Extension Cable StimConnect™		Surgical Pacing Electrode
	Connection Pacemaker- Side	Connection Patient-Side	
 <p>PACE® Model 203 with unipolar terminals accepting pins of 0.9 to 2.0 mm  (order number: P/N 2002)</p>	 <p><b>Article No: 81987BL (blue) / 81987 WS (white)</b></p> <p>Plugs with touch protected 2 mm pins fitting unipolar terminals of PACE® Model 203</p> <p>Connects exposed (naked) pins with diameter between 0.9 to 2.0 mm</p> <p>Packaging Unit: 5 Pieces</p>		 <p>Adapters with exposed (naked) 2 mm pins</p> <p>Exposed (naked) 0.9 mm pins/ breaking needles</p>
	 <p><b>Article No: 82207BL</b></p> <p>Plugs with exposed 2 mm pins fitting unipolar terminals of PACE® Model 203</p> <p>Covered alligator clips allow versatile connection</p> <p>Packaging Unit: 10 Pieces</p>		
	Reusable Extension Cable StimConnect™		
	 <p><b>Article No: 81828BL (blue) / 81828 WS (white)</b></p> <p>Plugs with exposed 2 mm pins fitting unipolar terminals of PACE® Model 203</p> <p>Connects exposed (naked) pins with diameter between 0 to 2.0 mm</p> <p>Packaging Unit: 5 Pieces</p>		 <p>Adapters with exposed (naked) 2 mm pins</p> <p>Exposed (naked) 0.9 mm pins/ breaking needles</p>

Please contact [mail@osypkamed.com](mailto:mail@osypkamed.com) for other StimConnect® extension cable variants.

**OSYPKA**  
MEDICAL

# StimConnect™ for PACE® Model 203 with Unipolar Terminals

## Extension Cables for Surgical Pacing Applications

PACE® Model 203 is primarily used for temporary surgical cardiac pacing applications. Surgical wires (temporary myocardial electrodes or heart wires) for temporary pacing are either equipped with 2 mm exposed (naked) plug pins or 0.9 mm pins/ needles. Following StimConnect® extension cables are recommended for PACE® Model 203 with bipolar terminals.

Disposable cables are for single use, reusable extension cables can be sterilized 25 times by ethylene oxide or steam. A punch card allows the identification of the sterilization cycles.

Pacemaker	Disposable Extension Cable StimConnect™		Surgical Pacing Electrode
	Connection Pacemaker- Side	Connection Patient-Side	
 <p>PACE® Model 203 with bipolar terminals (order number: P/N 2003)</p>	 <p><b>Article No: 82203BL (blue) / 82203WS (white)</b></p> <p>Bipolar plug fitting bipolar terminals of PACE® Model 203</p> <p>Connects exposed (naked) pins with diameter between 0.9 to 2.0 mm</p> <p>Packaging Unit: 10 Pieces</p>		 <p>Adapters with exposed (naked) 2 mm pins</p>  <p>Exposed (naked) 0.9 mm pins/breaking needles</p>
	 <p><b>Article No: 81832 WS</b></p> <p>Bipolar plug fitting bipolar terminals of PACE® Model 203</p> <p>Covered alligator clips allow versatile connection</p> <p>Packaging Unit: 10 Pieces</p>		
	Reusable Extension Cable StimConnect™		
	 <p><b>Article No: 82202BL (blue) / 82202WS (white)</b></p> <p>Bipolar plug fitting bipolar terminals of PACE® Model 203</p> <p>Connects exposed (naked) pins with diameter between 0 to 2.0 mm</p> <p>Packaging Unit: 5 Pieces</p>		 <p>Adapters with exposed (naked) 2 mm pins</p>  <p>Exposed (naked) 0.9 mm pins/breaking needles</p>

Please contact [mail@osypkamed.com](mailto:mail@osypkamed.com) for other StimConnect® extension cable variants.

 **OSYPKA**  
MEDICAL









# StimConnect™ for PACE® Model 300 with Unipolar Terminals

## Extension Cables for Surgical Pacing Applications

PACE® Model 300 is primarily used for temporary surgical cardiac pacing applications. Surgical wires (temporary myocardial electrodes or heart wires) for temporary pacing are either equipped with 2 mm exposed (naked) plug pins or 0.9 mm pins/ needles. Alligator clips facilitate extraordinary temporary pacing applications.

Disposable cables are for single use, reusable extension cables can be sterilized 25 times by ethylene oxide or steam. A punch card allows the identification of the sterilization cycles.

Pacemaker	Disposable Extension Cable StimConnect™		Surgical Pacing Electrode
	Connection Pacemaker- Side	Connection Patient-Side	
 <p>PACE® Model 300 with unipolar terminals accepting pins of 0.9 to 2.0 mm  (order number: P/N 2005)</p>	 <p><b>Article No: 81987BL (blue) / 81987 WS (white)</b></p> <p>Plugs with touch protected 2 mm pins fitting unipolar terminals of PACE® Model 203</p> <p>Connects exposed (naked) pins with diameter between 0.9 to 2.0 mm</p> <p>Packaging Unit: 5 Pieces</p>		 <p>Adapters with exposed (naked) 2 mm pins</p> <p>Exposed (naked) 0.9 mm pins/ breaking needles</p>
	 <p><b>Article No: 82207BL</b></p> <p>Plugs with exposed 2 mm pins fitting unipolar terminals of PACE® Model 203</p> <p>Covered alligator clips allow versatile connection</p> <p>Packaging Unit: 10 Pieces</p>		
	Reusable Extension Cable StimConnect™		
	 <p><b>Article No: 81828BL (blue) / 81828 WS (white)</b></p> <p>Plugs with exposed 2 mm pins fitting unipolar terminals of PACE® Model 203</p> <p>Connects exposed (naked) pins with diameter between 0 to 2.0 mm</p> <p>Packaging Unit: 5 Pieces</p>		 <p>Adapters with exposed (naked) 2 mm pins</p> <p>Exposed (naked) 0.9 mm pins/ breaking needles</p>

Please contact [mail@osypkamed.com](mailto:mail@osypkamed.com) for other StimConnect® extension cable variants.