


# EndoStim™

TEMPORARY PACING LEADS



# Temporary Pacing Leads

## EndoStim™

Pacing leads for transvenous applications to temporarily stimulate the ventricle or the atrium. Based on the application of external pacemakers suitable for emergency, pre- and postoperative stimulation as well as for rapid-pacing during TAVI procedures.

### Characteristics

The EndoStim™ portfolio includes bipolar and quadripolar pacing leads. A floating balloon catheter completes the product range.

The following French sizes, shapes, lengths, inter-electrode distances and connectors are available:

	Bipolar EndoStim™ Leads	Quadripolar EndoStim™ Leads	Pacing Catheter with Balloon
<b>French Size</b>	4F, 5F, 6F	4F, 5F, 6F	5F
<b>Length</b>	90 cm, 125 cm	90 cm, 125 cm	110 cm
<b>Shape</b>	Straight, Cournand, Josephson	Straight, Josephson	Cournand
<b>Inter-Electrode Distance</b>	2 mm, 5 mm, 10 mm	5 mm, 10 mm, 2-5-2 mm	10 mm
<b>Connector</b>	Exposed 2 mm pin with protective cap	Shrouded 2 mm pin	Exposed 2 mm pin with protective cap

### Shape

EndoStim™ pacing leads are available in three shapes: Straight, Cournand and Josephson.



### Connector

EndoStim™ pacing leads are either available with shrouded pins or with exposed pins with protective cap.



shrouded pins  
(Bipolar pacing leads)



exposed pins with protective cap  
(Quadripolar pacing leads, pacing lead with balloon)

## Order Information Temporary Pacing Leads

<b>Bipolar EndoStim™ Temporary Pacing Leads*</b>						
<b>(Packaging Unit: 5 pieces)</b>						
<b>Article No.</b>	<b>Number of Poles</b>	<b>French Size</b>	<b>Length</b>	<b>Shape</b>	<b>Inter-Electrode Distance</b>	<b>Connector</b>
27039-OMG	2	4F	90 cm	Straight	2 mm	Exposed Pins
27003-F	2	4F	90 cm	Straight	5 mm	Exposed Pins
27010-OMG	2	4F	90 cm	Straight	10 mm	Exposed Pins
27040-OMG	2	4F	90 cm	Cournand	2 mm	Exposed Pins
27018-OMG	2	4F	90 cm	Cournand	5 mm	Exposed Pins
27015-OMG	2	4F	90 cm	Cournand	10 mm	Exposed Pins
27038-OMG	2	4F	90 cm	Josephson	2 mm	Exposed Pins
27006-F	2	4F	90 cm	Josephson	5 mm	Exposed Pins
27009-OMG	2	4F	90 cm	Josephson	10 mm	Exposed Pins
27051-OMG	2	5F	90 cm	Straight	2 mm	Exposed Pins
27045-OMG	2	5F	90 cm	Straight	5 mm	Exposed Pins
27048-OMG	2	5F	90 cm	Straight	10 mm	Exposed Pins
27052-OMG	2	5F	90 cm	Cournand	2 mm	Exposed Pins
27046-OMG	2	5F	90 cm	Cournand	5 mm	Exposed Pins
27049-OMG	2	5F	90 cm	Cournand	10 mm	Exposed Pins
27050-OMG	2	5F	90 cm	Josephson	2 mm	Exposed Pins
27044-OMG	2	5F	90 cm	Josephson	5 mm	Exposed Pins
27047-OMG	2	5F	90 cm	Josephson	10 mm	Exposed Pins
27042-OMG	2	5F	125 cm	Straight	2 mm	Exposed Pins
27004-F	2	5F	125 cm	Straight	5 mm	Exposed Pins
27012-OMG	2	5F	125 cm	Straight	10 mm	Exposed Pins
27043-OMG	2	5F	125 cm	Cournand	2 mm	Exposed Pins
27019-OMG	2	5F	125 cm	Cournand	5 mm	Exposed Pins
27016-OMG	2	5F	125 cm	Cournand	10 mm	Exposed Pins
27041-OMG	2	5F	125 cm	Josephson	2 mm	Exposed Pins
27007-F	2	5F	125 cm	Josephson	5 mm	Exposed Pins
27011-OMG	2	5F	125 cm	Josephson	10 mm	Exposed Pins
27063-OMG	2	6F	90 cm	Straight	2 mm	Exposed Pins
27057-OMG	2	6F	90 cm	Straight	5 mm	Exposed Pins
27060-OMG	2	6F	90 cm	Straight	10 mm	Exposed Pins
27064-OMG	2	6F	90 cm	Cournand	2 mm	Exposed Pins
27058-OMG	2	6F	90 cm	Cournand	5 mm	Exposed Pins
27061-OMG	2	6F	90 cm	Cournand	10 mm	Exposed Pins
27062-OMG	2	6F	90 cm	Josephson	2 mm	Exposed Pins
27056-OMG	2	6F	90 cm	Josephson	5 mm	Exposed Pins
27059-OMG	2	6F	90 cm	Josephson	10 mm	Exposed Pins
27054-OMG	2	6F	125 cm	Straight	2 mm	Exposed Pins
27005-F	2	6F	125 cm	Straight	5 mm	Exposed Pins
27014-OMG	2	6F	125 cm	Straight	10 mm	Exposed Pins
27055-OMG	2	6F	125 cm	Cournand	2 mm	Exposed Pins
27020-OMG	2	6F	125 cm	Cournand	5 mm	Exposed Pins
27017-OMG	2	6F	125 cm	Cournand	10 mm	Exposed Pins
27053-OMG	2	6F	125 cm	Josephson	2 mm	Exposed Pins
27008-F	2	6F	125 cm	Josephson	5 mm	Exposed Pins
27013-OMG	2	6F	125 cm	Josephson	10 mm	Exposed Pins

\* Products and options may not be available or may vary by region and country. Please contact [mail@osypkamed.com](mailto:mail@osypkamed.com) for availability in your region and country and for more information.

# Temporary Pacing Leads

## Order Information Temporary Pacing Leads

Quadripolar EndoStim™ Temporary Pacing Leads*							(Packaging Unit: 5 pieces)
Article No.	Number of Poles	French Size	Length	Shape	Inter-Electrode Distance	Connector	
27022-OMG	4	4F	90 cm	Straight	5 mm	Shrouded Pins	
27028-OMG	4	4F	90 cm	Straight	10 mm	Shrouded Pins	
27033-OMG	4	4F	90 cm	Straight	2-5-2 mm	Shrouded Pins	
27021-OMG	4	4F	90 cm	Josephson	5 mm	Shrouded Pins	
27027-OMG	4	4F	90 cm	Josephson	10 mm	Shrouded Pins	
27032-OMG	4	4F	90 cm	Josephson	2-5-2 mm	Shrouded Pins	
27024-OMG	4	5F	125 cm	Straight	5 mm	Shrouded Pins	
27030-OMG	4	5F	125 cm	Straight	10 mm	Shrouded Pins	
27035-OMG	4	5F	125 cm	Straight	2-5-2 mm	Shrouded Pins	
27023-OMG	4	5F	125 cm	Josephson	5 mm	Shrouded Pins	
27029-OMG	4	5F	125 cm	Josephson	10 mm	Shrouded Pins	
27034-OMG	4	5F	125 cm	Josephson	2-5-2 mm	Shrouded Pins	
27026-OMG	4	6F	125 cm	Straight	5 mm	Shrouded Pins	
27065-OMG	4	6F	125 cm	Straight	10 mm	Shrouded Pins	
27037-OMG	4	6F	125 cm	Straight	2-5-2 mm	Shrouded Pins	
27025-OMG	4	6F	125 cm	Josephson	5 mm	Shrouded Pins	
27031-OMG	4	6F	125 cm	Josephson	10 mm	Shrouded Pins	
27036-OMG	4	6F	125 cm	Josephson	2-5-2 mm	Shrouded Pins	

Bipolar EndoStim™ Temporary Pacing Lead with Balloon*								(Packaging Unit: 1 piece)
Article No.	Number of Poles	French Size	Length	Shape	Inter-Electrode Distance	Ballon Diameter	Ballon Volume	Connector
29207-BS	2	5F	110 cm	Cournand	10 mm	9 mm	1.0 cc	Exposed Pins



Latex balloon

## EasyStim™ Adapter – Accessory for Temporary Pacing Leads

The EasyStim™ adapter allows to connect shrouded 2mm pins to collets for exposed (nacked) 2mm pins.

EasyStim™ Adapter*			(Packaging Unit: 20 double packs)
Article No.	Product Name	Sterility	
94005	EasyStim	no-sterile	
94006	EasyStim.S	sterile	



\* Products and options may not be available or may vary by region and country. Please contact [mail@osypkamed.com](mailto:mail@osypkamed.com) for availability in your region and country and for more information.

Osypka Medical GmbH

Albert-Einstein-Strasse 3  
 12489 Berlin, Germany  
 Phone: +49 (30) 6392 8300  
 Fax: +49 (30) 6392 8301  
 E-Mail: [mail@osypkamed.com](mailto:mail@osypkamed.com)  
[www.osypkamed.com](http://www.osypkamed.com)



EasyStim and EndoStim are trademarks of Osypka Medical.  
 © Osypka Medical 2024