






StimConnect™

for PACE® Model 101 with Unipolar Terminals

Extension Cables for Transvenous Pacing Applications

PACE® Model 101 is primarily used for temporary transvenous cardiac pacing applications. Transvenous leads for temporary pacing are either equipped with 2 mm shrouded plug pins or exposed (naked) plug pins. Following StimConnect® extension cables are recommended for PACE® Model 101 with unipolar terminals.

Pacemaker	Disposable Extension Cable StimConnect™		Transvenous Pacing Lead
	Connection Pacemaker- Side	Connection Patient-Side	
 <p>PACE® Model 101 with unipolar terminals accepting pins of 0.9 to 2.0 mm (order number: P/N 2001).</p>	 <p>Article No: 82207BL</p> <p>Touch protected plug fitting unipolar terminals of PACE® Model 101</p> <p>Connects shrouded pins of transvenous leads</p> <p>Packaging Unit: 10 Pieces</p>	 <p>Shrouded pins (ø 2 mm)</p>	
	 <p>Article No: 81987BL</p> <p>Touch protected plug fitting unipolar terminals of PACE® Model 101</p> <p>Connects exposed (naked) pins of transvenous leads</p> <p>Packaging Unit: 5 Pieces</p>	 <p>Exposed (naked) pins (ø 2 mm)</p>	

OSYPKA
MEDICAL

Products and options may not be available or may vary by region and country. Please contact mail@osypkamed.com for availability in your region and country, for more information and for other StimConnect® extension cable variants.

StimConnect™

for PACE® Model 101 with Bipolar Terminal

Extension Cables for Transvenous Pacing Applications

PACE® Model 101 is primarily used for temporary transvenous cardiac pacing applications. Transvenous leads for temporary pacing are either equipped with 2 mm shrouded plug pins or exposed (naked) plug pins. Following StimConnect® extension cables are recommended for PACE® Model 101 with bipolar terminals.

Pacemaker	Disposable Extension Cable StimConnect™		Transvenous Pacing Lead
	Connection Pacemaker- Side	Connection Patient-Side	
 <p>PACE® Model 101 with bipolar terminal (order number: P/N 2004)</p>	 <p>Article No: 82208BL</p> <p>Bipolar plug fitting bipolar terminal of PACE® Model 101</p> <p>Packaging Unit: 10 Pieces</p>	 <p>Article No: 82203BL</p> <p>Bipolar plug fitting bipolar terminal of PACE® Model 101</p> <p>Packaging Unit: 10 Pieces</p>	 <p>Shrouded pins (ø 2 mm)</p>
	 <p>Exposed (naked) pins (ø 2 mm)</p>		

OSYPKA
MEDICAL

Products and options may not be available or may vary by region and country. Please contact mail@osypkamed.com for availability in your region and country, for more information and for other StimConnect® extension cable variants.